

Nucleus® Troubleshooting Kit for Educators

Cochlear Americas is pleased to offer the Nucleus® Troubleshooting Kit for Educators. Designed to assist those who work with children with Nucleus cochlear implants in school settings, this kit was developed to help identify and remedy simple problems with a child's sound processor and minimize the time a child is without auditory stimulation.

The troubleshooting kit itself is a canvas bag that holds spare components for a specific processing system. There are kits available for each of the following Nucleus sound processors:

- Nucleus® 22 Spectra Bodyworn
- Nucleus® 22 ESPrin™ 22 ear-level
- Nucleus® 22 ESPrin™ 3G
- Nucleus® 24 SPrin™ Bodyworn
- Nucleus® 24 ESPrin™ ear-level
- Nucleus® 24 ESPrin™ 3G ear-level
- Nucleus® Freedom™ BTE
- Nucleus® Freedom™ Bodyworn

Each kit comes with a manual which takes the reader through a step-by-step process to troubleshoot in various areas. The manual, which is sectioned according to each sound processing system, includes a detailed illustration of the system as well as a list of kit contents, instructions for performing a maintenance check, and step-by-step troubleshooting tips for common problem areas specific to that system. Troubleshooting tips include how to change batteries and how to check the transmitting coil, headset, and microphone. The guide also features helpful illustrations and diagrams.



Each Nucleus® Troubleshooting Kit is ordered based upon the external sound processor the child is utilizing. Pricing information is as follows. Please call Cochlear Customer Service at 1-877-833-3101 to place your order.

**FUZI52 Troubleshooting Kit \$ 180.00
for SPrint™ sound processor**

Includes: 4" HS8 transmitting cable, 24" HS8 head set cable, large earhook, 31" lapel microphone & clip, signal check, monitor earphones, 2 AA batteries and a Troubleshooting Guide, all in a nylon zipper case.

**FUZI53 Troubleshooting Kit \$215.00
for ESPrit™ sound processor,
for Nucleus® 24 recipients**

Includes: 5" ESPrit transmitting cable, 675 Power One zinc air batteries (6 pack), signal check, large earhook, 59" Lapel microphone & clip, monitor earphones and a Troubleshooting Guide, all in a nylon zipper case.

**FUZI54 Troubleshooting Kit \$195.00
for ESPrit™ 3G sound processor**

Includes: 675 Power One zinc air batteries (6 pack), soft earhook cover, earmold compatible earhook, 59" Lapel microphone & clip, monitor earphones, signal check, accessories adapter and a Troubleshooting Guide, all in a nylon zipper case.

**FUZI55 Troubleshooting Kit \$140.00
for Spectra sound processor**

Includes: 24" HS6 Headset Cable, 4" Transmitting Cable, earhooks for techmin mic, 31" lapel microphone, signal check, 2 AA batteries and a Troubleshooting Guide, all in a nylon zipper case.

**FUZI56 Troubleshooting Kit \$215.00
for ESPrit™ 22 sound processor**

Includes: 4" ESPrit 22 transmitting cable, signal check, monitor earphones, earhook, 675 Power One zinc air batteries (6 pack), 59" lapel microphone clip and a Troubleshooting Guide, all in a nylon zipper case.

**FUZI70 Troubleshooting Kit \$315.00
for Freedom™ BTE sound processor**

Includes: 675 Power One zinc air batteries (6 pack), Freedom coil (8cm), magnet X2, lapel microphone, monitor earphones, earhook set, microphone protector, battery holder and a Troubleshooting Guide, all in a nylon zipper case.

**FUZI71 Troubleshooting Kit \$375.00
for Freedom™ Bodyworn
sound processor**

Includes: Energizer disposable AAA batteries (2), Freedom coil (8cm), magnet X2, lapel microphone, monitor earphones, earhook set, microphone protector, bodyworn shoe and cable and a Troubleshooting Guide, all in a nylon zipper case.

The Nucleus Troubleshooting Kit for Educators is an effective support tool for clinicians and educators who work with children with Nucleus cochlear implants, offering valuable information that can save time and spare the silence that could result from a minor problem.

www.cochlear.com

Cochlear Americas
400 Inverness Parkway
Suite 400
Englewood CO 80112
USA

Cochlear Limited
14 Mars Road
Lane Cove NSW 2066
Australia

Tel: 1 303 790 9010
Fax: 1 303 792 9025
Toll Free: 1 800 523 5798

Tel: 61 2 9428 6555
Fax: 61 2 9428 6352
Toll Free: 1 800 620 929

Nucleus is a registered trademark of Cochlear Limited. Freedom, SPrint, and ESPrit are trademarks of Cochlear Limited.

FUZI79 ISSUE 1 AUG05